



Provation Apex Q3 2023 Release Highlights

Persistent Pathology Results Filters

- Filters applied to the Pathology Results grid will now be retained when loading in new result sets.

Expanded text limit in Data Reports

- The character limit for displaying results in Data Reports has been increased to 300 characters.

Interface Message History: Show full HL7 Messages

- Full HL7 messages will be displayed when viewing previously sent messages.

Patient Merge

- In the Manage Patients settings, users can now identify possible duplicates and merge patient records together.

Audit Log

- A new report was added to the Provation CSV reports. The Audit Log report will allow users to track significant apex events over a specific period.

Amend any Procedure Note

- Under the Exam Access tab, specific security groups can be permitted to put finalized notes into amendment mode 60 days after the original procedure date.

Delete Exams

- Under the Exam Access tab, specific security groups can be granted the ability to Delete Exams, including patient information and documentation.

Rename Groups in Security Settings

- Users with select access can click on the new pencil icon in the title bar and rename the existing security groups.

Image Capture Additional Security

- A new security setting was added to allow select users access to Image Capture without giving them full access to the procedure note.

Date of Birth and Order Number added to Patient Profile Search

- Apex has added the ability to search by patient Date of Birth or Order Number.