



Provation® Apex

## Offline Image Capture

Your safety net – ensuring productivity, even if you lose connectivity.

Intelligent cloud-based Provation® Apex relies on internet connection to allow clinicians to document procedures, embed captured images and attest to findings.

When network connectivity is lost before or during a procedure, even temporarily, clinicians still have the ability to save images locally on their computer using Provation Apex's Offline Image Capture feature.

Offline Image Capture is built to ensure that loss of network communication does not interfere with patient care. If network issues occur, clinicians can enter Offline Mode, select their current patient's ID and continue taking images.

Offline Image Capture helps you:

- ✓ Mitigate the risk of lost procedure images
- ✓ Minimize the need for delayed exams
- ✓ Ensure a positive patient experience
- ✓ Reduce reliance on an organization's network connection

When connection is restored, images can be automatically associated to the correct patient procedure in Provation Apex.

### With Offline Image Capture you can:

- ✓ Capture all relevant images
- ✓ Use free text to create temporary patient entries
- ✓ Access protected health information (PHI)\*

\*If the clinician has provided their authenticated credentials while online before Offline Image Capture is launched.

