



Provation® iPro HL7 Interfaces and Connections with MEDITECH

Maximize your IT investment and improve workflow efficiency

The Provation iPro platform leverages HL7 interfaces to integrate with and extend existing IT system investments. Provation has made it easy for facilities to streamline workflows and support the seamless exchange of health information.

Benefits:


- Maximizes utilization of existing IT systems
- Enhance clinician productivity by eliminating redundancy from clinical workflows
- Get experienced, professional implementation
- Support quality of patient care and patient experience

INBOUND

From MEDITECH Expance to Provation iPro

Scheduling Information	SIU
Patient Demographics	ADT
Allergies	ADT
Preop vitals	ORU or PCS
Preop labs	ORU
Auto vitals	ORU
Home Medications	RAS
Inpatient Medications	RDE
Preop Health Issues	FHIR

Provation® iPro Platform + Meditech Expance



OUTBOUND

From Provation iPro to MEDITECH Expance

MDM or SFTP Cold Feed

PDF of Completed Anesthesia Record from iPro

DFT	Anesthesia Billing Data
ADT	Allergies
FHIR	Anesthesia Administered Medications
RAS	Home Medications
PPR	Preop Health Issues

